

Anti-IFNGR1 Antibody (A11871)


Specifications:

Name:	Anti-IFNGR1 Antibody
Description:	Rabbit polyclonal antibody to IFNGR1.
Applications:	WB, IHC
Recommended Dilutions:	WB: 1:500-1:1,000, IHC: 1:50-1:200
Reactivity:	Human, Mouse, Rat
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 270-489 of human IFNGR1 (NP_000407.1).
Sequence:	IKKINPLKEKSIILPKSLISVVRSATLETKPESKYVSLITSYQPFSLEKEVVCEEPLS PATVPGMHTEDNPGKVEHTEELSSITEVVTTEENIPDVVPGSHLTPIERESSSPLSSN QSEPGSIALNSYHSRNCSESDHSRNGFDTDSSCLESHSSLSDSEFPPNNKGEIKTEGQ ELITVIKAPTSFGYDKPHVLVDLLVDDSGKESLIGYRPTEDSKEFS
l la atr	
Host:	Rabbit
Host: Clonality:	Rabbit Polyclonal
Clonality:	Polyclonal
Clonality: Isotype:	Polyclonal IgG
Clonality: Isotype: Conjugate:	Polyclonal IgG Unconjugated
Clonality: Isotype: Conjugate: Purification:	Polyclonal IgG Unconjugated Affinity purification.
Clonality: Isotype: Conjugate: Purification: Molecular Weight:	Polyclonal IgG Unconjugated Affinity purification. 80 - 90 kDa
Clonality: Isotype: Conjugate: Purification: Molecular Weight: Product Form:	Polyclonal IgG Unconjugated Affinity purification. 80 - 90 kDa Liquid


antibodies

Anti-IFNGR1 Antibody (A11871)


Images:



Western blot analysis of extracts of various cell lines, using Anti-IFNGR1 Antibody (A11871) at 1:1,000 dilution. The secondary antibody was Goat Anti-Rabbit IgG H&L Antibody (HRP) at 1:10,000 dilution. Lysates/proteins were present at 25µg per lane. The blocking buffer used was 3% non-fat dry milk in TBST. Detection was with a ECL Basic Kit. Exposure time: 180s.



Immunohistochemistry analysis of paraffin-embedded mouse spleen using Anti-IFNGR1 Antibody (A11871) at a dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of paraffin-embedded rat spleen using Anti-IFNGR1 Antibody (A11871) at a dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.