

# **Anti-MyD88 Antibody (A89417)**

### Specifications:

Name: Anti-MyD88 Antibody

Description: Rabbit polyclonal antibody to MyD88.

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1:500-1:1,000, ICC/IF: 1:50-1:200

Reactivity: Human, Mouse, Rat

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human

MyD88 (NP\_002459.3).

Sequence: LAEEMDFEYLEIRQLETQADPTGRLLDAWQGRPGASVGRLLELLTKLGRDDVLLELGP

SIEEDCQKYILKQQQEEAEKPLQVAAVDSSVPRTAELAGITTL

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Conjugate: Unconjugated

Purification: Affinity purification.

Molecular Weight: 33 kDa


Product Form: Liquid

Formulation: Supplied in Phosphate Buffered Saline, pH 7.3, with 50% Glycerol and 0.01% Thiomersal.

Storage: Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

Disclaimer: This product is for research use only. It is not intended for diagnostic or therapeutic use.

#### Images:




Western blot analysis of extracts of Mouse spleen, using Anti-MyD88 Antibody (A89417) at 1:1,000 dilution. The secondary antibody was Goat Anti-Rabbit IgG H&L Antibody (HRP) at 1:10,000 dilution. Lysates/proteins were present at 25µg per lane. The blocking buffer used was 3% non-fat dry milk in TBST. Detection was with a ECL Basic Kit. Exposure time: 1s.




# Anti-MyD88 Antibody (A89417)


## Images continued:



Western blot analysis of extracts of Rat spleen, using Anti-MyD88 Antibody (A89417) at 1:1,000 dilution. The secondary antibody was Goat Anti-Rabbit IgG H&L Antibody (HRP) at 1:10,000 dilution. Lysates/proteins were present at 25µg per lane. The blocking buffer used was 3% non-fat dry milk in TBST. Detection was with a ECL Basic Kit. Exposure time: 60s.



Immunofluorescence analysis of A549 cells using Anti-MyD88 Antibody (A89417) at a dilution of 1:200 (40x lens). DAPI was used to stain the cell nuclei (blue).



Immunofluorescence analysis of C6 cells using Anti-MyD88 Antibody (A89417) at a dilution of 1:200 (40x lens). DAPI was used to stain the cell nuclei (blue).