

Anti-RPL36 Antibody (A15900)

Specifications:

Name: Anti-RPL36 Antibody

Description: Rabbit polyclonal antibody to RPL36.

Applications: WB, IHC, ICC/IF

Recommended Dilutions: WB: 1:500-1:2,000, IHC: 1:50-1:200, ICC/IF: 1:50-1:200

Reactivity: Human, Mouse, Rat

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of

human RPL36 (NP_056229.2).

Sequence: MALRYPMAVGLNKGHKVTKNVSKPRHSRRGRLTKHTKFVRDMIREVCGFAPYERRAM

ELLKVSKDKRALKFIKKRVGTHIRAKRKREELSNVLAAMRKAAAKKD

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Conjugate: Unconjugated

Purification: Affinity purification.

Molecular Weight: 15 kDa

Product Form: Liquid


Formulation: Supplied in Phosphate Buffered Saline, pH 7.3, with 50% Glycerol and 0.02% Sodium

Azide.

Storage: Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

Disclaimer: This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:




Western blot analysis of extracts of various cell lines, using Anti-RPL36 Antibody (A15900) at 1:1,000 dilution. The secondary antibody was Goat Anti-Rabbit IgG H&L Antibody (HRP) at 1:10,000 dilution. Lysates/proteins were present at 25µg per lane. The blocking buffer used was 3% non-fat dry milk in TBST. Detection was with a ECL Basic Kit. Exposure time: 10s.




Anti-RPL36 Antibody (A15900)


Images continued:



Immunohistochemistry analysis of paraffin-embedded rat brain using Anti-RPL36 Antibody (A15900) at a dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of paraffin-embedded human liver cancer using Anti-RPL36 Antibody (A15900) at a dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of paraffin-embedded mouse heart using Anti-RPL36 Antibody (A15900) at a dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.