

Anti-MLX Antibody [PCRP-MLX-1G8] (A250082)

Specifications:

Name: Anti-MLX Antibody [PCRP-MLX-1G8]

Description: Mouse monoclonal [PCRP-MLX-1G8] antibody to MLX.

Applications: Flow Cytometry, IF, WB

Recommended Dilutions: Flow Cytometry: 1-2 μg/million cells, IF: 1-2 μg/ml, WB: 1-2 μg/ml

Reactivity: Human

Immunogen: Recombinant full-length human MLX protein.

Host: Mouse

Clonality: Monoclonal

Clone ID: PCRP-MLX-1G8

Isotype: IgG2b

Conjugate: Unconjugated

Purification: Protein A/G chromatography.

Concentration: 200 µg/ml

Product Form: Liquid

Formulation: Supplied in 10mM Phosphate Buffered Saline with 0.05% BSA and 0.05% Sodium Azide.

Storage: Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

General Notes: This monoclonal antibody is also available in a different formulation without BSA and


Sodium Azide - Anti-MLX Antibody [PCRP-MLX-1G8] - BSA and Azide free (A253262).

Disclaimer: This product is for research use only. It is not intended for diagnostic or therapeutic use.




Anti-MLX Antibody [PCRP-MLX-1G8] (A250082)


Images:



Flow cytometric analysis of PFA fixed HeLa cells using Anti-MLX Antibody [PCRP-MLX-1G8] followed by Goat Anti-Mouse IgG (CF® 488) (Blue). Unstained cells (red).



Immunofluorescent analysis of PFA fixed HeLa cells stained with Anti-MLX Antibody [PCRP-MLX-1G8] followed by Goat Anti-Mouse IgG (CF® 488) (Green). Counterstain is Phalloidin (Red).




SDS-PAGE analysis of Anti-MLX Antibody [PCRP-MLX-1G8] under non-reduced and reduced conditions; showing intact IgG and intact heavy and light chains, respectively. SDS-PAGE analysis confirms the integrity and purity of the antibody.



Anti-MLX Antibody [PCRP-MLX-1G8] (A250082)

Images continued:



Analysis of protein array containing more than 19,000 full-length human proteins using Anti-MLX Antibody [PCRP-MLX-1G8]. Z-Score and S- Score: The Z-score represents the strength of a signal that a monoclonal antibody (MAb) (in combination with a fluorescently-tagged anti-IgG secondary antibody) produces when binding to a particular protein on the HuProtTM array. Z-scores are described in units of standard deviations (SD's) above the mean value of all signals generated on that array. If targets on HuProtTM are arranged in descending order of the Z-score, the S-score is the difference (also in units of SD's) between the Z-score. S-score therefore represents the relative target specificity of a MAb to its intended target; a MAb is considered to be specific to its intended target, if the MAb has an S-score of at least 2.5. For example, if a MAb binds to protein X with a Z-score of 43 and to protein Y with a Z-score of 14, then the S-score for the binding of that MAb to protein X is equal to 29.